

## PAIN RELIEF AND ALIGNMENT

The patented Hallufix (R) Hallux valgus splint was developed by scientists and orthopedic surgeons for pain relief and alignment of the widely common malpositioned 1st toe (Hallux).

## PREVENTIVE AND POSTOPERATIVE

Preventively, mild and medium severe Hallux valgus deformities benefit from the Hallufix splint. A deterioration of the Hallux' misalignment can be delayed or even stopped. After a Hallux valgus surgery the splint helps to secure the result of the operation and to shorten the rehabilitation time.

## EASY TO USE

The Hallufix® splint is easy to apply and can be worn left or right. As of the hinge joint, the device can be worn as either a night splint or even during the day under load (when walking). Wearing the splint in shoes with wide forefoot area is presumed – UXGO shoes are the perfect match. The incorporated gel pad provides cushion and can be cooled to ease the pain of the bunion.



## AT A GLANCE

### APPLICATION / PURPOSE


Hallux valgus deformity

### BENEFITS

- Counteracts the malaligned Hallux at rest, when walking and standing
- Prevents deterioration of the deformity
- Can reduce postsurgical rehabilitation time
- Maintains foot mobility
- Protects and cools the Hallux' sensitive joint area
- Tailored to the outstanding UXGO footwear from Perpedes

### ORDER INFORMATION

ITEM-NO.	DESCRIPTION
HALCS	HALLUFIX, Comfort splint
HALCSXL	HALLUFIX, Comfort splint, XL
HALNS	HALLUFIX, Replacement kit



Order and support:  
 Mon - Thu: 8.00 am - 5.30 pm  
 Fri: 8.00 am - 4.00 pm  
 T. +49 (0)7021 738 30-0  
 F. +49 (0)7021 738 30-298  
 info@perpedes.de

Return address:  
 Perpedes GmbH  
 Tannenbergr. 139  
 73230 Kirchheim/Teck  
 Germany

2018/09/01



**Hallufix®**  
 INNOVATIVE ORTHOTIC SOLUTIONS

**HALLUX VALGUS  
 SPLINT**

CONSERVATIVE ALIGNMENT  
 INSTEAD OF SURGERY

# STAY MOBILE WITH THE HALLUFIX SPLINT



## A SUITABLE SHOE FOR SENSITIVE FEET

Worldwide countless women suffer from a Hallux valgus. The search for a suitable shoe with sufficient volume in the forefoot area for painless walking is tiring. In cooperation with Hallufix, PERPEDES developed the outstanding UXGO footwear that is not merely „wide fit“ but provides adequate support for the rearfoot combined with a roomier forefoot area. This allows a pain-free, natural gait pattern with or without a Hallufix splint.

UXGO sneakers are insole-friendly and available in two widths – „COMFORT FIT“ and „COMFORT PLUS FIT“.

**CHOOSE FROM MORE THAN 30 STYLES AND FEEL THE DIFFERENCE!**




Browse online in the UXGO catalog.




## OPERATION PRINCIPLE

WITHOUT HALLUFIX



$\alpha$  Hallux valgus angle  
 $\beta$  Intermetatarsal angle

WITH HALLUFIX



Hallufix splint

### MIDFOOT SPLINT

Connects the brace to the midfoot. Incorporated gel pad for cooling purposes



### TOE SPLINT

Re-aligns the Hallux' malposition



### INTEGRATED HINGE JOINT

Allows forefoot mobility

Gel pad